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 Achtung!  **LOT** Chargennummer  **REF** Artikelnummer

Instructions For Use Pelvic Base Training Aid Phase II

Dear customer,

Thank you for choosing our ELANEE pelvic base training aid – phase II! These instructions will help you utilize this product as best as possible. Please read the instructions for use thoroughly before the first training with the ELANEE pelvic base training aid – phase II!

We are at your disposal for any questions or suggestions you may have and wish you success with your pelvic base training!

Your ELANEE team

 Please observe instructions for use and packaging and keep them!

Made in Germany

1. Field of application

These instructions for use apply to the following product:

REF 705-00 ELANEE pelvic base training aid phase II

2. Product content

- 1 Training aids, 1 Instructions for use

3. Purpose and intended use

The pelvic base training aid – phase II supports the **permanent follow-up training**, after you have already learned to perceive your pelvic base and a certain core strength of the muscles has been built (e.g. with the ELANEE pelvic base training aids – phase I).

- **During a pregnancy without problems**, you can train with the pelvic base training aid – phase II for as long as you are feeling comfortable with it.
- **After childbirth**, you can begin with the training as soon as

the cervix has closed and your pelvic muscles are strengthened again. In case of doubt, please contact your midwife or your physician.

Pelvic base training with the ELANEE pelvic base training aid – phase II is therefore suitable:

- For permanent training of the pelvic base muscles with an already pre-trained pelvic base (e.g. with the ELANEE pelvic base training aids – phase I).
- For increasing sexual sensations.

4. Functioning

The ELANEE pelvic base training aid – phase II is composed of two tubular bodies that are interconnected by a wide strip. Both tubular bodies contain a ball floating freely. The internal balls swing against the outside wall of the tubular body with every motion and generate slight vibration. This vibration is directly transferred to the pelvic base muscles. This stimulates and activates the muscles.

The large supporting surfaces of the lower part help transfer the vibrations directly to the pelvic base muscles. The 14 imprinted stimulation islands facilitate additional stimulus to the muscles. The wide strip between both halves enhances the swinging of both balls and the training effect is intensified by this.

The ELANEE pelvic base training aid – phase II is completely inserted. To simplify the insertion, the top part of the cone is designed smaller than the bottom part. It is invisible from the outside and silent. The training ensues by simply wearing the training aid and can therefore be incorporated in everyday life. The pelvic base muscles are stimulated by every motion.

5. Application

1. Wash the ELANEE pelvic base training aid – phase II thoroughly before the first use and make sure that conditions are clean (e.g. wash hands).
2. Insert the small part of the cone in the front first. The application of lube glide gel is recommended if the vagina is dry when training with the pelvic base training aid – phase II.
3. Insert the training aid slowly, completely and approximately as far as you would insert a tampon. Both crescent-shaped supporting surfaces should be right and left during insertion. This way, they rest directly on the pelvic base muscles.
4. Try in the beginning how long you feel wearing them is comfortable and increase the length according to your well-being.
5. Conduct the training daily so that your pelvic base muscles remain active.

6. Safety information

-  Thoroughly clean the ELANEE Pelvic Floor Training Aids both prior to initial usage and after every application and ensure clean conditions (e.g. wash hands).
-  Please check your training aids before initial usage. Use only your cones if they are in perfect condition. If the training aid is scratched or otherwise damaged this can cause injury.
-  In cases of vaginal infections, inflammations, injury or painful complaints in the genital area and during menstruation or a problematic pregnancy, you should avoid using the Pelvic Floor Training Aids.

-  When other complaints (such as a dropping of the pelvic organs) exist, you should have other causes than a weakness of the pelvic floor ruled out in advance by a doctor before you start with the training.
-  In cases of a known hypersensitivity and allergies to one of the named ingredients, you should not use the aids or only use them after consultation with a doctor.
-  If you are a virgin, you should not use the pelvic floor training aids.
-  The wearing of cervical caps, contraceptive rings, diaphragms or intrauterine devices during training are possible. However, you must check the correct positioning before and after each training session. If in doubt only use after consultation with a doctor.
-  Do not use steam sterilisation, dry heat sterilisation or boiling to clean the aids as this can lead to distortion of the material.

7. Cleaning information

- Clean the training aid with warm water and customary soap or disinfect the cone with a customary skin disinfection agent.
- Before using the training aid again, possible residue of soap or disinfectant should be rinsed off thoroughly using water.
- Gas sterilisation is also possible.

8. Disposal information

- You can dispose of the ELANEE pelvic base training aid phase II in the household waste.

- Please separate the individual components of the packaging correctly sorted and recycle them in the appropriate recycling system.
- Please observe local or regional disposal regulations.

9. Material composition

The ELANEE pelvic base training aids phase II consist of:

- Polypropylene
- Thermoplastic elastomer
- 2 rubberised steel balls

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